
Cricothyrotomy Trainers

Directions for Use

Deluxe Model #135



 **Simulaids**

Assembly:

This unit comes ready to use, no assembly is required.

Deluxe Unit

Replacing Overlay Skin:

1. Remove airway.
2. Pull skin off base.
3. Replace new skin over base.
4. Replace airway.

Replacing the Airway:

1. Remove overlay skin (see above).
2. Remove trachea.
3. Pull airway out of base.
4. Reinsert (child or adult) airway into base. (Make sure airway is fully inserted in base.)
5. Replace matching trachea.
6. Pull skin taut, reattach velcro (on hand held model).

Cleaning:

Use any mild detergent to clean outside of trainer.

Replacement Parts:

#137 Adult Trachea

#138 Child Trachea

#132 Adult Airway (4)

#133 Child Airway (4)

#131 Handheld Overlay Skins (2)

#136 Deluxe Overlay Skins (4)

#PP1366 Carry Bag

Return Policy

Should it be necessary to return an item for any reason, contact our Customer Service Department to obtain an RGA Number. Please refer to your invoice number when phoning in your request for returning merchandise. Should you have any questions or wish further information on any product we manufacture call or write our Customer Service Department.

Warranty

Simulaids warrants this product to be free from any defect in materials and/or workmanship for a period of three years from the date of purchase, as evidenced by the date of invoice when the product was shipped to the end user. This warranty expressly does not cover abuse, accidental or purposeful damage, or any form of modification to the product. Simulaids reserves the right to either repair or replace affected parts or the entire unit, at their sole discretion, after investigating and reviewing the actual product and the damage. In most instances, a digital photo of the product in question showing the damage will help qualify a product for return to the factory. At no time will any product be accepted at the plant without proper return authorization issued by Simulaids.

Freight and Shipping charges are the sole responsibility of the end user. No product will be received with shipping charges due. Any product considered for warranty work must be identified by serial number and invoice number from the agency through whom the product was purchased. Without this information the product will not receive a return authorization number as required above.



 **Simulaids**



PO Box 1289

16 Simulaids Dri

Saugerties NY 12477

Toll Free: (800)-431-4310 Fax: (845)-679-8996

www.simulaids.com

3/2006